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INSTRUCTIONS FOR CUSTOM MADE SHOULDER ARTICULAR PROSTHESES

GENERAL PRODUCT INFORMATION

Even though the surgeon can now count on the means necessary to restore mobility, correct deformity, and reduce pain for many patients, thanks to the advancements in partial and total joint replacement, which generally allow him to achieve these goals successfully, it must be remembered that the prostheses used are manufactured from metal, ceramic and plastic materials and that any joint replacement system, therefore, cannot be expected to withstand activity levels and loads to avoid normal healthy tone. In addition, the system will not have the same strength, reliability or duration as a natural human joint.

- In using total joint prostheses, the surgeon should be aware of the following:
 - The correct selection of the prosthesis is extremely important.
 - The potential for success in total joint replacement is increased by selection of the proper size, shape, and design of the prosthesis. Total joint prostheses require careful seating and adequate bone support. Smaller sized implants are intended for patients with small bone and normally slight weight. Such components could be inappropriate for other patients. Surgeons are encouraged to use their best medical judgement when choosing the proper implant size regardless of the endosteal area of the bone.
 - In selecting patients for total joint replacements, the following factors can be critical to the eventual success of the procedure.
 - Patient's weight. An overweight or obese patient can produce high loads on the prosthesis, which can lead to failure of the prosthesis. This becomes a major consideration when the patient is small boned and a small size prosthesis must be used.
 - Patient's occupation or activity. If the patient is involved in an occupation or activity, which calls for substantial walking, running, lifting, or muscle strain, the resultant forces can lead to failure of the fixation. The device or both. The prosthesis will not restore function to the level expected with normal healthy bone, and the patient should not have unrealistic functional expectations.
 - Condition of sensory, mental, illness, or alcoholism. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the prosthesis, leading to failure or other complications.
 - Foreign body sensitivity. In case of suspected material sensitivity, appropriate tests should be made prior to material selection or implantation.

Caution: Patients with renal insufficiency may be sensitive to potential metal ion release. Also, since little is known about the transport of potential metal ion release across the placenta, these devices should be used with caution in women of childbearing age.

GENERAL

The biocompatible materials used in the manufacturing of implantable medical devices show an excellent resistance to corrosion and conform to specific international standards:

- Stainless steel, standard ISO 5832 - 1 or 9, symbol: Stainless steel
- Titanium alloy, standard ISO 5832 - 3, symbol: Ti6Al4V
- Cobalt-Chrome-Molybdenum alloy, standard ISO 5832 - 4 or 12, symbol: Co-Cr-Mo

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- Notes
- * Symbol ISO 7000: Graphical Symbols For Use On Equipment - Registered Symbols
- 2. From symbol #1 to symbol #13 in accordance with ISO 15223-1
- 3. Symbol #15 in accordance with MDD 93/42/EEC

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- Alumina ceramics, standard ISO 6474, symbol: Al2O3
- Ultra high molecular weight polyethylene, standard ISO 5834 - 1 and 2, symbol: UHMWPE
- Unalloyed titanium, standard ISO 5832 - 2, symbol: Ti. The symbol of the material used to manufacture the device is specified on the product label.

IMPLANTATION WITH OTHER MEDICAL DEVICES

The various components of implantable medical devices manufactured by Adler Ortho are designed in such a way as to ensure their interchangeability and compatibility. In particular:

- If ceramic heads are intended to be used, modular necks and mono-block stems manufactured by Adler Ortho should only be employed in combination with ceramic heads manufactured by Adler Ortho. The Company performs specific tests on these heads, to guarantee a safe coupling to the Adler Ortho taper.
- Metal / Metal inserts have to be used only in combination with Metal / Metal heads manufactured by Adler Ortho.

Adler Ortho take no responsibility for the association of their implants with components from another manufacturer.

INAPPROPRIATE COMBINATIONS OF MATERIALS IN ARTICULATING SURFACES

- Stainless steels / titanium alloys
- Stainless steels / stainless steels
- Stainless steels /unalloyed titanium
- Stainless steel (ISO 5832 - 1) / cobalt chrome alloys
- Unalloyed titanium / unalloyed titanium
- Unalloyed titanium / titanium alloys
- Unalloyed titanium / cobalt chrome alloys
- Unalloyed titanium / ultra high molecular weight polyethylene
- Titanium alloys / cobalt chrome alloys
- Ceramics/metal
- Zirconia ceramics /alumina ceramics

INAPPROPRIATE COMBINATIONS OF METALS IN NON ARTICULATING CONTACT SURFACES

- Stainless steel (excluding the stainless steel described in ISO 5832 - 9) / cobalt chrome alloys
- Stainless steel (excluding the stainless steel described in ISO 5832 - 9) / unalloyed titanium

INDICATIONS

- Primary and secondary arthrosis
- Cuff tear arthropathy
- Rheumatoid and degenerative arthritis
- Deformities, such as dysplasia
- Fractures and necrosis of bone
- Revisions, where other devices or treatments have failed
- Metastatic diseases

CONTRA-INDICATIONS

- Contra-indications may be relative or absolute. The surgical problems must be evaluated case by case, taking into account alternative surgical options (osteotomy, arthrodesis, amputation, etc.). The following examples are regarded as contra-indications:

- infection, septicæmia, and osteomyelitis constitute cases of absolute contra-

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- indicatie;
- serious metabolic, cardiovascular, respiratory or neurological pathologies;
- serious osteoporosis;
- rapid disease progression as manifested by joint destruction or bone absorption apparent on roentgenogram;
- obesity;
- skeletally immature patients;
- female patients of childbearing age, for whom a negative pregnancy test is not obtained;
- high patient activity which could lead to overloading of the implant.

NOTE

Given their complexity, it is advised that, in the selection of an artificial implant, the surgeon studies the clinical indications and contra-indications of that particular implant and that the surgeon uses his personal experience and judgement.

POSSIBLE COMPLICATIONS, SIDE-EFFCTS

Apart from possible per-operation complications, a prosthetic implant can also be subjected to:

- Peri-prosthetic infection with or without loosening
- Movement of one or several prosthetic elements due to mechanical overloading, osteoporosis, etc.
- Dislocation of the prosthesis or fracture of the bone due to trauma from physiotherapy
- Extra-articular pathology: phlebotrombosis, pulmonary embolism, etc.
- fatigue fracture of prosthetic components (e.g. as a result of : trauma, strenuous activity, improper alignment or duration of service.

VIGILANCE SYSTEM (EUROPEAN DIRECTIVE 93/42/EEC)

The users are party to the setting up of the vigilance system after introduction of medical device on the market. They are obliged by law to report to the competent authorities as well as the manufacturer any incident concerning the medical device liable to cause or having caused serious risk for the patient or user health.

VERIFICATION PRIOR TO USE - PRECAUTIONS

It is essential to have all the sizes listed in the catalogue in the operating theatre prior to the surgery to ensure the optimum size for the patient is selected. Presentation cases and instruments are available to achieve this. Maintenance and handling must be carried out in accordance with the recommendations of standard ISO 8828.

Handling must be limited and carried out with extreme care (non-powered gloves must be worn for the manipulation of devices coated with hydroxyapatite) so as to avoid damage/scratching, shock or contact with materials liable to alter the surface of the implant). In all cases the final conformation of the condition of the implant lies with the surgeon. Special care must be taken with ceramic devices which are not to be re-used even when in the absence of any apparent damage.

Only use a polyacetal tip to introduce ceramic devices.

At implantation all foreign bodies such as pieces of tissue and bone or cement particles must be removed. They could cause abnormal wear of the articulating surfaces.

If Morse tapers are used to assemble the implants, confirmation of the compatibility of the taper angles must be made. Before the neck is inserted, the female cone of the stem should be thoroughly cleaned. All components should be thoroughly cleaned before their assembly.

Use the femoral prostheses for dysplasia only in special cases, taking into account that the body weight may affect the mechanical resistance of the prosthesis

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