ADLER® ORTHO

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ADLER ORTHO

INSTRUCTIONS FOR CUSTOM MADE EN KNEE ARTICULAR PROSTHESES

GENERAL PRODUCT INFORMATION

Even though the surgeon can now count on the means necessary to restore mobility, correct deformity, and reduce pain for many patients, thanks to the advancements in partial and total joint replacement, which generally allow him to achieve these goals successfully, it must be remembered that the prostheses used are manufactured from metal and/or plastic materials and that any joint replacement system, therefore, cannot be expected to withstand activity levels and loads as would normal healthy bone. In addition, the system will not have the same strength, reliability or duration as a natural human joint.

- In using total joint prostheses, the surgeon should be aware of the following:
 - The correct selection of the prosthesis is extremely important.
 - The potential for success in total joint replacement is increased by selection of the proper size, shape, and design of the prosthesis. Total joint prostheses require careful seating and adequate bone support. Smaller sized implants are intended for patients with small bone and normally slight weight. Such components could be inappropriate for other patients. Surgeons are encouraged to use their clinical evaluation when choosing the proper implant size regardless of the endosteal area of the bone.
 - In selecting patients for total joint replacements, the following factors can be critical to the eventual success of the procedure.
1. Patient's weight. An overweight or obese patient can produce high loads on the prosthesis which can lead to failure of the prosthesis. This becomes a major consideration when the patient is small boned and a small size prosthesis must be used.
2. Patient's occupation or activity. If the patient is involved in an occupation or activity, which calls for substantial walking, running, lifting, or muscle strain, the resultant forces can lead to failure of the fixation, the device or both. The prosthesis will not restore function to the level expected with normal healthy bone, and the patient should not have unrealistic expectations.
3. Condition of senility, mental illness, or alcoholism. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the prosthesis, leading to failure or other complications.
4. Foreign body sensitivity. In case of suspected material sensitivity, appropriate tests should be made prior to material selection or implantation.

GENERAL

The biocompatible materials used in the manufacturing of implantable medical devices show an excellent resistance to corrosion and conform to specific international standards:

- Stainless steel, standard ISO 5832 - 1 or 9, symbol: Stainless steel
- Titanium alloy, standard ISO 5832 - 3, symbol: Ti6Al4V

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- Notes
1. * Symbol ISO 7000: Graphical Symbols For Use On Equipment - Registered Symbols
 2. From symbol #1 to symbol #13 in accordance with ISO 15223-1
 3. Symbol #15 in accordance with MDD 93/42/EEC

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- Cobalt-Chrome-Molibdenum alloy, standard ISO 5832 - 4 or 12, symbol: Co-Cr-Mo
- Alumina ceramics, standard ISO 6474, symbol: Al2O3
- Ultra high molecular weight polyethylene, standard ISO 5834 - 1 and 2, symbol: UHMWPE
- Unalloyed titanium, standard ISO 5832 - 2, symbol: Ti. The symbol of the material used to manufacture the device is specified on the product label.

IMPLANTATION WITH OTHER MEDICAL DEVICES
The various components of implantable medical devices manufactured by Adler Ortho are designed in such a way as to ensure their interchangeability and compatibility.

Adler Ortho take no responsibility for the association of their implants with components from another manufacturer.

- INAPPROPRIATE COMBINATIONS OF MATERIALS IN ARTICULATING SURFACES**
- Stainless steels / titanium alloys
 - Stainless steels / stainless steels
 - Stainless steels /unalloyed titanium
 - Stainless steel (ISO 5832 - 1) / cobalt chrome alloys
 - Unalloyed titanium / unalloyed titanium
 - Unalloyed titanium / titanium alloys
 - Unalloyed titanium / cobalt chrome alloys
 - Unalloyed titanium / ultra high molecular weight polyethylene
 - Titanium alloys / cobalt chrome alloys
 - Ceramics/metal
 - Zirconia ceramics /alumina ceramics

- INAPPROPRIATE COMBINATIONS OF METALS IN NON ARTICULATING CONTACT SURFACES**
- Stainless steel (excluding the stainless steel described in ISO 5832 - 9) / cobalt chrome alloys
 - Stainless steel (excluding the stainless steel described in ISO 5832 - 9) / unalloyed titanium

- INDICATIONS**
- Primary and secondary arthrosis
 - Rheumatoid arthritis
 - Bone necrosis
 - Revisions, where other devices or treatments have failed

CONTRA-INDICATIONS
Contra-indications may be relative or absolute. The articular problems must be evaluated case by case, taking into account alternative surgical options (osteotomy, arthrodesis, amputation, etc.). The following elements are regarded as contra-indications:

- infection, septicaemia, and osteomyelitis constitute cases of absolute contra-indication;
- serious metabolic, cardiovascular, respiratory or neurological pathologies;
- serious osteoporosis;
- rapid disease progression as manifested by joint destruction or bone absorption

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- ion apparent on roentgenogram;
- obesitas;
- skeletaly immature patients;
- femate patients of childbearing age, for whom a negative pregnancy test is not obtained;
- high patient activity which could lead to overloading of the implant.

NOTE
Given their complexity, it is advised that, in the selection of an articular implant, the surgeon studies the clinical indications and contra-indications of that particular implant and that the surgeon uses his personal experience and judgement.

POSSIBLE COMPLICATIONS, SIDE-EFFECTS
Apart from possible per-operation complications, a prosthetic implant can also be subject to:

- Peri-prosthetic infection with or without loosening
- Movement of one or several prosthetic elements due to mechanical overloading, osteoporosis, etc.
- Dislocation of the prosthesis or fracture of the bone due to traumatism from physiotherapy
- Extra-articular pathology: phlebotrombosis, pulmonary embolism, etc.
- Fatigue fracture of prosthetic components can occur as a result of : trauma, strenuous activity, improper alignment or duration of service.

VIGILANCE SYSTEM (EUROPEAN DIRECTIVE 93/42/EEC)
The users are party to the setting up of the vigilance system after introduction of medical device on the market. They are obliged by law to report to the competent authorities as well as the manufacturer any incident concerning the medical device liable to cause or having caused serious risk for the patient or user health.

VERIFICATION PRIOR TO USE - PRECAUTIONS
It is essential to have all the sizes listed in the catalogue in the operating theatre prior to the surgery to ensure the optimum size for the patient is selected. Presentation cases and instruments are available to achieve this. Maintenance and handling must be carried out in accordance with the recommendations of standard ISO 8828. Handling must be limited and carried out with extreme care (non-powered gloves must be worn for the manipulation of devices coated with hydroxyapatite) so as to avoid damage/scratching, shock or contact with materials liable to alter the surface of the implant). In all cases the final confirmation of the condition of the implant lies with the surgeon.

At implantation all foreign bodies such as pieces of tissue and bone or cement particles must be removed. They could cause abnormal wear of the articulating surfaces.

- If these papers are used to assemble the implants, confirmation of the dimensional compatibility of the taper must be made.
- All components should be thoroughly cleaned before their assembly. Under no circumstances should a damaged or scratched prosthesis or a prosthesis having undergone any inappropriate or unauthorised treatments be implanted, this includes implants that have already been used.
- Never steam sterilize/resterilize plastic components.
- Please check the sterilization by verifying the colour of the sterilization in-

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indicator (red in case of gamma sterilization, green in case of ethylene oxide sterilization). Please verify the vacuum tightness of the packaging of the knee prosthesis metal components. Do not use a device in which the packaging has been damaged or is questionable because it could call into question the integrity of the packaging and the sterility of the device. Please verify the sterilization expiration date.

RE-USE
Implants must never be reused as their mechanical and/or functional characteristics could be definitively compromised so their service life could be reduced.

INSTRUMENTS SUPPLIED WITH THE PROSTHESIS
Instruments supplied with the prosthesis must be sterilized by means of a validated and controlled sterilization method.

For metallic instruments we recommend steam sterilization (Autoclave) usually employed in hospitals. 134°C cycles for 18 minutes do not cause any deformation or damage to the plastic or metal instrument.

It is important to avoid any contact with objects or corrosive substances (e.g.: strong acid solutions, bases, solvents) which could adversely alter the material.

DEVICES SUPPLIED IN A STERILE STATE
The metallic components of the knee prosthesis supplied by Adler Ortho in a sterile state are sterilized in vacuum by means of gamma rays radiation at a minimum dose of 25 KGy. The polyethylene components of the knee prosthesis supplied by Adler Ortho in a sterile state are sterilized by means of ethylene oxide. The label clearly shows the sterilisation indicator and the indication "STERILE".

RESTERILISATION
All implants supplied by Adler Ortho must not be re-sterilized by the purchaser. Adler Ortho does not take any responsibilities for the use of implants re-sterilized by the purchaser.

STORAGE CONDITIONS
The implants must be stored under dry conditions and be protected from sunlight, extremes in temperature and contamination by particles.

PATIENT INFORMATION
The patient destined to receive the implant must be informed that the security and durability of the implant depends also on his/her behaviour, particularly as concerning his/her activity and body weight.

For any further information please contact your supplier.

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