# ADLER ORTHO

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## INSTRUCTIONS FOR CUSTOM MADE ELBOW PROSTHESIS

### GENERAL PRODUCT INFORMATION

Even though the surgeon can now count on the means necessary to restore mobility, correct deformity, and reduce pain for many patients, thanks to the advancements in partial and total joint replacement, which generally allow him to achieve these goals successfully, it must be remembered that the prostheses used are manufactured from metal, ceramic and plastic materials and that any joint replacement system, therefore, cannot be expected to withstand activity levels and loads as would normal healthy bone. In addition, the system will not have the same strength, reliability or duration as a natural human joint.

In using elbow prosthesis, the surgeon should be aware of the following:

- The correct selection of the prosthesis is extremely important.
- The potential for success is increased by selection of the proper size, shape, and design of the prosthesis.
- Joint prostheses require careful seating and adequate bone support.

When choosing the proper implant, surgeons are encouraged to use their best medical judgement.

In selecting patients for elbow arthroplasty, the following factors can be critical to the eventual success of the procedure.

1. Patient's occupation or activity. If the patient is involved in an occupation which calls for a high level of activity, weightlifting, or muscle strain, the resultant forces can lead to failure of the fixation, the device or both. The prosthesis will not restore function to the level expected with normal healthy bone, and the patient should not have unrealistic functional expectations.

2. Condition of senility, mental illness, or alcoholism. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the prosthesis, leading to failure or other complications.

3. Foreign body sensitivity. In case of suspected material sensitivity, appropriate tests should be made prior to material selection or implantation.

**Caution:** Patients with renal insufficiency may be sensitive to potential metal ion release. Also, since little is known about the transport of potential metal ion release across the placenta, these devices should be used with caution in women of childbearing age.

### GENERAL

The biocompatible materials used in the manufacturing of implantable medical devices show an excellent resistance to corrosion and conform to specific international standards:

- Stainless steel, standard ISO 5832 - 1 or 9, symbol: Stainless steel
- Titanium alloy, standard ISO 5832 - 3, symbol: Ti6Al4V
- Cobalt-Chrome-Molibdenum alloy, standard ISO 5832 – 4 or 12, symbol: Co-Cr-Mo

-Alumina ceramic, standard ISO 6474, simbolo: Al2O3

-Ultra high molecular weight polyethylene, standard ISO 5834 – 1 and 2, simbolo: UHMWPE

-Unalloyed titanium, standard ISO 5832 - 2, simbolo: Ti

The symbol of the material used to manufacture the device is specified on the

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product label.

### IMPLANTATION WITH OTHER MEDICAL DEVICES

The various components of implantable medical devices manufactured by Adler Ortho are designed in such a way as to ensure their interchangeability and compatibility.

Adler Ortho take no responsibility for the association of their implants with components from another manufacturer.

### INAPPROPRIATE COMBINATIONS OF MATERIALS IN ARTICULATING SURFACES

- Stainless steels / titanium alloys
- Stainless steels / stainless steels
- Stainless steels /unalloyed titanium
- Stainless steel (ISO 5832 - 1) / cobalt chrome alloys
- Unalloyed titanium / unalloyed titanium
- Unalloyed titanium / titanium alloys
- Unalloyed titanium / cobalt chrome alloys
- Unalloyed titanium / ultra high molecular weight polyethylene
- Titanium alloys / cobalt chrome alloys
- Ceramics/metal
- Zirconia ceramic /alumina ceramics

### INAPPROPRIATE COMBINATIONS OF METALS IN NON ARTICULATING CONTACT SURFACES

- Stainless steel (excluding the stainless steel described in ISO 5832 - 9) / cobalt chrome alloys
- Stainless steel (excluding the stainless steel described in ISO 5832 - 9) / unalloyed titanium

### INDICATIONS

- Fracture type III or IV of Mason's classification
- Elbow fracture, luxation
- Degenerative joint diseases
- Post-traumatic bone necrosis

### CONTRA-INDICATIONS

Contra-indications may be relative or absolute.

The articular problems must be evaluated case by case, taking into account alternative surgical options (osteosynthesis, excision of radial epiphysis, amputation, et.).

The following examples are regarded as contra-indications:

- infection, septicæmia, and osteomyelitis constitute cases of absolute contra-indication;
- serious metabolic, cardiovascular, respiratory or neurological pathologies;
- serious osteoporosis;
- rapid disease progression as manifested by joint destruction or bone absorp2tion apparent on roentgenogram;
- neurogenic arthropathy (Charcot joint);
- skeletally immature patients;
- femala patients of childbearing age, for whom a negative pregnancy test is not

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