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en INSTRUCTIONS FOR HIP ARTICULAR PROSTHESES - RECTA STEM

Even though the surgeon can now count on the means necessary to restore mobility, correct deformity, and reduce pain for many patients, thanks to the advancements in partial and total joint replacement, which generally allow him to achieve these goals successfully, it must be remembered that the prostheses used are manufactured from metal, ceramic and plastic materials and that any joint replacement system, therefore, cannot be expected to withstand activity levels and loads as would normal healthy bone. In addition, the system will not have the same strength, reliability or duration as a natural human joint.

- In using total joint prostheses, the surgeon should be aware of the following:
- The correct selection of the prosthesis is extremely important.
 - The potential for success in total joint replacement is increased by selection of the proper size, shape, and design of the prosthesis. Total joint prostheses require careful seating and adequate bone support. Smaller sized implants are intended for patients with small bone and normally slight weight. Such components could be inappropriate for other patients. Surgeons are encouraged to use their best medical judgement when choosing the proper implant size regardless of the endosteal area of the bone.
 - In selecting patients for total joint replacements, the following factors can be critical to the eventual success of the procedure.

1. Patient's weight. An overweight or obese patient can produce high loads on the prosthesis, which can lead to failure of the prosthesis. This becomes a major consideration when the patient is small boned and a small size prosthesis must be used.
2. Patient's occupation or activity. If the patient is involved in an occupation or activity, which calls for substantial walking, running, lifting, or muscle strain, the resultant forces can lead to failure of the fixation, the device or both. The prosthesis will not restore function to the level expected with normal healthy bone, and the patient should not have unrealistic functional expectations.
3. Condition of senility, mental illness, or alcoholism. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the prosthesis, leading to failure or other complications.
4. Foreign body sensitivity. In case of suspected material sensitivity, appropriate tests should be made prior to material selection or implantation. Caution: Patients with renal insufficiency may be sensitive to potential metal ion release. Also, since little is known about the transport of potential metal ion release across the placenta, these devices should be used with caution in women of childbearing age.

GENERAL

The biocompatible materials used in the manufacturing of implantable medical devices show an excellent resistance to corrosion and conform to specific international standards:

- Stainless steel, standard ISO 5832 - 1 or 9, symbol: Stainless steel
 - Titanium alloy, standard ISO 5832 - 3, symbol: Ti6Al4V
 - Cobalt-Chrome-Molybdenum alloy, standard ISO 5832 - 4 or 12, symbol: Co-Cr-Mo
 - Alumina ceramics, standard ISO 6474, symbol: Al2O3
 - Ultra high molecular weight polyethylene, standard ISO 5834 - 1 and 2, symbol: UHMWPE
 - Unalloyed titanium, standard ISO 5832 - 2, symbol: Ti
- The symbol of the material used to manufacture the device is specified on the product label.

IMPLANTATION WITH OTHER MEDICAL DEVICES

The various components of implantable medical devices manufactured by Adler Ortho are designed in such a way as to ensure their interchangeability and compatibility.

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In particular:

- if ceramic heads are intended to be used, modular necks and mono-block stems manufactured by Adler Ortho should only be employed in combination with ceramic heads manufactured by Adler Ortho. The Company performs specific tests on these heads, to guarantee a safe coupling to the Adler Ortho taper.
- Metal / Metal inserts have to be used only in combination with Metal / Metal heads manufactured by Adler Ortho.

The RECTA femoral stem can only be used with the following types of modular necks manufactured by Adler Ortho: 0X, 0A, 0Y, 9X, 9Y, 9A, 9AA

Adler Ortho take no responsibility for the association of their implants with components from another manufacturer.

INAPPROPRIATE COMBINATIONS OF MATERIALS IN ARTICULATING SURFACES

- Stainless steels / Titanium alloys
- Stainless steels / stainless steels
- Stainless steels / stainless steels
- Stainless steels /un alloyed titanium
- Stainless steel (ISO 5832 - 1) / cobalt chrome alloys
- Unalloyed titanium / unalloyed titanium
- Unalloyed titanium / titanium alloys
- Unalloyed titanium / cobalt chrome alloys
- Unalloyed titanium / ultra high molecular weight polyethylene
- Titanium alloys / cobalt chrome alloys
- Ceramics/metal
- Zirconia ceramics /alumina ceramics

INAPPROPRIATE COMBINATIONS OF METALS IN NON ARTICULATING CONTACT SURFACES

- Stainless steel (excluding the stainless steel described in ISO 5832 - 9) / cobalt chrome alloys
- Stainless steel (excluding the stainless steel described in ISO 5832 - 9) / unalloyed titanium

INDICATIONS

- Primary and secondary arthritis
- Rheumatoid and degenerative arthritis
- Deformities, such as dysplasia
- Fractures and necrosis of bone
- Revisions, where other devices or treatments have failed
- Metastatic diseases

CONTRA-INDICATIONS

- Contra-indications may be relative or absolute.
- The articular problems should be evaluated case by case, taking into account alternative surgical options (osteotomy, arthrodesis, amputation, etc.)
- The following examples are regarded as contra-indications:
- infection, septicemia, and osteomyelitis constitute cases of absolute contra-indication;
 - serious metabolic, cardiovascular, respiratory or neurological pathologies;
 - serious osteoporosis;
 - rapid disease progression as manifested by joint destruction or bone absorption apparent on roentgenogram;
 - obesity;
 - skeletal maturation patients;
 - female patients of childbearing age, for whom a negative pregnancy test is not

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