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**en** STERILE SINGLE-USE SURGICAL INSTRUMENTS

**PRODUCT GENERAL INFORMATION**

These use instructions apply to all sterile single-use surgical instruments manufactured by Adler Ortho. Using these instruments requires the knowledge of biomechanical anatomy and reconstructive surgery of joints, thus they must be used solely by a qualified surgeon with the necessary expertise and experience in the field of joint replacement and solely for the purposes described in the surgical technique.

**GENERAL FEATURES**

The biocompatible material used to manufacture the sterile instruments is the Nylon. The indication of the material making up the medical device appears on the label applied to the packaging.

**STERILE DEVICES**

The surgical instruments are delivered in sterile conditions by Adler Ortho. They are vacuum sterilized by gamma ray irradiation with a minimum dose of 25 kGy. A red toning stamp and the word "STERILE" clearly appear on the label.

**CHECK BEFORE OPERATION. CAUTIONS**

Check the conditions of the instruments before use. Check sterilization by making sure that the toning stamp is of red color. Check that the packaging is vacuum-packed.

Perform handling carefully and thoroughly to avoid any deterioration. Do not use the instrument whose packaging is damaged or of poor quality as the packaging integrity and therefore the sterility of the device might be compromised. Check the expiry date of the sterilization before use.

The use of the instruments can generate heat caused by the friction between instruments and bones. It is recommended to cool down the instruments by using physiological saline solution. After use, instruments must be carefully inspected to check the absence of any damage and to prevent any residue from entering the patient's body. After resection, wash the bone with physiological solution to avoid any residue in the wound.

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- Notes  
1. \* Symbol ISO 7000: Graphical Symbols For Use On Equipment - Registered Symbols  
2. From symbol #1 to symbol #13 in accordance with ISO 15223-1  
3. Symbol #15 in accordance with MDD 93/42/EEC

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Instruments or parts thereof must not be implanted under any circumstances.

**STORAGE AND HANDLING**

Store the devices in a dry room away from sunlight, extreme temperatures and dust. Incorrect handling can damage the instruments. When handling the packaged instruments, care must be taken not to damage the packaging.

**DISPOSAL**

Sterile single-use instruments must be treated and disposed of as common single-use medical devices contaminated with biological material.

**WARNINGS**

Sterile single-use instruments must never be reused as their mechanical and/or functional characteristics might be permanently compromised and their useful life reduced. The purchaser must never sterilize the sterile single-use instruments. Adler Ortho does not assume any responsibility in the event of devices reesterilized by the purchaser.

**MONITORING (EUROPEAN DIRECTIVE 93/42/EEC)**

Users are concerned party in the application of the monitoring system after placing the medical device on the market. They are required by law to notify the competent authorities and the manufacturers of any incidents relating to the medical device which have entailed or might entail a serious risk for the health of the patient or a user.

For any additional information contact the supplier.

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